

## Creatine Kinase (CK) Diagnostic Reagent Kit (DGKC Method) Product Instruction

### Product name

Creatine Kinase (CK) Diagnostic Reagent Kit  
(DGKC Method)

Abbreviation: CK

### Package Specifications

| R1     | R2     |
|--------|--------|
| 80ml×2 | 20ml×1 |

Calibrator (lyophilized):

Specification 1 (1.0mL×1; 1 level).

### Intended use


Creatine kinase (CK) diagnostic reagent is used for quantitative in vitro determination of creatine kinase in human serum on photometric systems.

### Clinical significance

The determination of CK and CK-isoenzyme activities is utilized in the diagnosis and monitoring of myocardial infarction such as the progressive Duchenne muscular dystrophy. Following injury to the myocardium, such as occurs with acute myocardial infarction, CK is released from the damaged myocardial cells. In early cases, a rise in the CK-activity can be found just 4 hour after an infarction. The CK-activity reaches a maximum after 12-24 hours and then falls back to the normal range after 3-4 days.

### Principle

The assay consists of three reaction steps:



The NADPH formed is measured by determining its absorbance at 340 nm and this is directly proportional to the amount of CK activity in the sample.

### Main ingredient

| Composition        | Concentration |
|--------------------|---------------|
| Imidazole buffer   | 100mmol/L     |
| N-acetyl cysteine  | 20mmol/L      |
| EDTA               | 2mmol/L       |
| NADP               | 2mmol/L       |
| G6P-DH             | >1500U/L      |
| Hexokinase (HK)    | >2500U/L      |
| Glucose            | 20mmol/L      |
| Creatine phosphate | 30mmol/L      |

|                                 |          |
|---------------------------------|----------|
| Di(adenosine-5') pentaphosphate | 10μmol/L |
| ADP                             | 2mmol/L  |

### Storage and stability

This kit stays stable for 12 months when stored in a dark condition at 2-8°C (do not freeze). After being opened for the first time, the kit can be stored at 2-8°C in an anti-pollution environment and stays stable for 14 days.

### Specimen preparation

Serum.

Stability:

2 days at 15-25°C

7 days at 2-8°C

4 weeks at (-15)-(-25)°C.

### Assay procedure

|                                                                                                                                       | Blank  | Sample |
|---------------------------------------------------------------------------------------------------------------------------------------|--------|--------|
| Reagent 1                                                                                                                             | 200μl  | 200μl  |
| Dist. water                                                                                                                           | 10μl   | -      |
| Sample                                                                                                                                | -      | 10μl   |
| Mix. Incubate for 5 minutes at 37°C. Then add:                                                                                        |        |        |
| Reagent 2                                                                                                                             | 50μl   | 50μl   |
| Mix. Incubate for 2 minutes at 37°C. Read absorbance within 2 minutes and calculate the average rate of the absorbance change ΔA/min. |        |        |
| Main wavelength                                                                                                                       | 340 nm |        |
| Sub wavelength                                                                                                                        | 405 nm |        |
| Method                                                                                                                                | Rate   |        |
| Reaction direction                                                                                                                    | +      |        |

### Procedure (Manual):

1. Zero spectrophotometer at 340nm with distilled water.
  2. Appropriately label cuvettes: reagent blank, sample, etc.
  3. Pipette 1.0ml (1000μl) of reagent 1 into all cuvettes.
  4. Add 0.50ml (50μl) of sample to respective cuvettes. Mix.
  5. Incubate for 5 minutes at 37°C.
  6. Then add 0.25ml (250μl) of reagent 2 to respective cuvettes. Mix gently.
  7. Incubate for 2 minutes at 37°C. Read absorbance within 2 minutes and calculate the average rate of the absorbance change ΔA/min.
- Determine the average absorbance per minute (ΔA/min), multiply by factor 4180 for results in U/L.

### Calculations

$$\text{Conc. of CK} = \frac{A(\text{sample}) - A(\text{blank})}{A(\text{Cal.}) - A(\text{blank})} \times \text{conc. of Cal.}$$

### Reference interval

Males: < 190 U/L

Females: < 170 U/L

It is recommended that each laboratory should establish its own reference interval.

### Quality control

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### Performance characteristics

Measuring range

The method is linear up to 1000 U/L.

Samples above this concentration should be diluted 1+1 with 0.9% NaCl solution and the result multiplied by 2.

### Interference

The effect of the following substances can be neglected if the concentrations of the following substances are at or below the given values.

| Substances       | Concentrations |
|------------------|----------------|
| Direct Bilirubin | 40mg/dl        |
| Total Bilirubin  | 40mg/dl        |
| Haemoglobin      | 200mg/dl       |
| Intralipid       | 0.50 %         |
| VC               | 50mg/dl        |

### Automation

Special adaptations for automatic analyzers can be made on request.

### General precautions

For in vitro diagnostic use only.

Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines.

### Precautions for measurement

Specimens should be treated as potentially infectious (HIV, Hepatitis B virus, Hepatitis C virus, etc.) and handled with appropriate caution. Reagents with different lot numbers should not be interchanged or mixed.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Material safety data sheet available for professional user on request.

### Parameters

| HITACHI Parameters |                      |
|--------------------|----------------------|
| Model              | 7180                 |
| Assay Code         | Rate-A               |
| Reaction Time      | 10                   |
| Assay Point        | 22-32                |
| Wave Length Sec    | 405                  |
| Wave Length Prim   | 340                  |
| ABS. Limit         | 32000                |
| Sample Volume (uL) | 10                   |
| R1 Volume (uL)     | 200                  |
| R2 Volume (uL)     | 0                    |
| R3 Volume (uL)     | 50                   |
| Calibration Type   | Linear               |
| Calibration Point  | 2                    |
| Span Point         | 2                    |
| SD Limit           | 999                  |
| Duplicate Limit    | 1000 or (0%-1000Abs) |
| Sensitivity Limit  | 0                    |
| Prozone Limit      | 32000                |

| OLYMPUS 400        |           |
|--------------------|-----------|
| Sample Volume (uL) | 10        |
| Sample DIL. VOL.   | /         |
| R1 Volume (uL)     | 200       |
| R2 Volume (uL)     | 50        |
| Reagent DIL. VOL.  | /         |
| Measure Wave 1     | 340       |
| Measure Wave 2     | 405       |
| METHOD             | RATE      |
| SLOPE              | -         |
| Method Point       | 15--25    |
| Linear Limit       | -         |
| NO LAG TIME        | Yes       |
| OD L               | -         |
| OD H               | -         |
| Reagent OD Limit   |           |
| FIRST              | -2.0-+2.5 |

|                      |           |
|----------------------|-----------|
| LAST                 | -2.0-+2.5 |
| Dynamic Range L      | 15        |
| Dynamic Range H      | 600       |
| Correlation Factor A | 1.0       |
| Correlation Factor B | 0.0       |
| Calibration Type     | AB        |
| Formula              | liner     |
| STD(1) Conc.#        | *         |
| Unit                 | U/L       |

**Dubai - Global Office****Xportin Global FZE**

Business Centre,  
Sharjah Publishing City Free Zone,  
Sharjah, United Arab Emirates  
Formation No. 4307793.01  
Mobile: +971 506938802 / +91 98241 88802  
Email: [xportinglobalfzelle@gmail.com](mailto:xportinglobalfzelle@gmail.com)  
Website: [www.xportin.com](http://www.xportin.com)

**India - Head Office:****Devasvi Enterprise LLP (Manufacturer)**

Navsari, Gujarat – 396445  
Email: [devashvienterprisellp@gmail.com](mailto:devashvienterprisellp@gmail.com)