

## Immunoglobulin G (IgG) Diagnostic Reagent Kit (Immunoturbidimetry Method) Product Instruction

### Product name

Generic Name: Immunoglobulin G (IgG)  
Diagnostic Reagent Kit (Immunoturbidimetry Method)  
Abbreviation: IgG

### Package Specifications

| R1         | R2         |
|------------|------------|
| 30ml × 1   | 10ml × 1   |
| 750ml × 1  | 250ml × 1  |
| 7500ml × 1 | 2500ml × 1 |

Calibrator (lyophilized):  
Specification 1 (1.0mL × 4; 4 levels);  
Quality Control (lyophilized):  
Specification 1 (1.0mL × 1; 1 level).

### Intended use

Immunoglobulin G (IgG) diagnostic reagent is used for the quantitative in vitro determination of IgG in human serum on photometric systems.

### Clinical significance

IgG is the main component of immunoglobulin in serum, accounting for about 75% of the total content of immunoglobulin in serum. Synthesized mainly by plasma cells in the spleen and lymph nodes, it is the only antibody that can cross the placenta and plays a major role in preventing infection in the first few weeks of life in newborns. The content of IgG varies greatly among individuals, and it fluctuates also greatly by the same individual under different conditions. Most of the antibacterial, antiviral, and antitoxin antibodies produced by the body under the stimulation of antigens belong to IgG, which is the main antibody of the body's anti-infection immunity. Commonly used detection methods are immunoturbidimetry and chemiluminescence.

### Principle

When Immunoglobulin G in human serum and its corresponding antibody (goat anti-human immunoglobulin G antibody) meet in the liquid phase, an antigen-antibody complex immediately forms and a certain turbidity also forms. Compared with the calibrators treated in the same way, the content of immunoglobulin G in the sample can be calculated.

### Main ingredient

| Composition | Main ingredients             |
|-------------|------------------------------|
| R1          | Phosphate buffer 20mmol/L    |
|             | Lipid elimination agent 5g/L |
|             | Sodium chloride 9g/L         |

|                 |                                            |
|-----------------|--------------------------------------------|
| R2              | Goat anti-human IgG specific antiserum 33% |
| Calibrator      | Tris 50mM                                  |
|                 | Trehalose 5%                               |
|                 | BSA 5%                                     |
|                 | Human immunoglobulin G                     |
| Quality Control | Tris 50mM                                  |
|                 | Trehalose 5%                               |
|                 | BSA 5%                                     |
|                 | Human immunoglobulin G                     |

### Storage and stability

This kit stays stable for 18 months when stored in a dark condition at 2-8°C (do not freeze). After being opened for the first time, the kit can be stored at 2-8°C in an anti-pollution environment and stays stable for 14 days.

### Specimen preparation

Fresh serum;  
-20°C, stable for 7 days.

### Assay procedure

#### 1. Reagent preparation

This reagent is liquid and can be used directly.

|                      |       |               |               |
|----------------------|-------|---------------|---------------|
| Main wavelength      | 570nm | Sample(S)     | 3 μl          |
| Subwavelength        | 800nm | Reagent 1     | 225 μl        |
| Reaction temperature | 37°C  | Reagent 2     | 75 μl         |
| Cuvette light path   | 1cm   | Reaction type | Two Point End |

#### 2. Operation steps

|                                                                                                                                                |        |
|------------------------------------------------------------------------------------------------------------------------------------------------|--------|
| Add to the cuvette:                                                                                                                            |        |
| Sample (S)                                                                                                                                     | 3 μl   |
| R1                                                                                                                                             | 225 μl |
| Mix well, incubate at 37°C for 5 minutes, 570nm, set zero with a blank tube, and measure the absorbance A1, then add                           |        |
| R2                                                                                                                                             | 75 μl  |
| Mix well, incubate at 37°C for 5 minutes, 570nm, set zero with a blank tube, and measure the absorbance A2, absorbance A2 $\Delta A = A2 - A1$ |        |

### Calculations

Calculate and plot  $\Delta A = (A2 - A1)$  of the calibrators versus assigned concentration values on a linear-linear graph paper. Calculate  $\Delta A$  optical densities of samples and read values in mg/dl on the reference curve. Samples yielding absorbances above highest calibrator should be retested after further dilution.

**Reference interval**

8.00 g/L -16.00g/L

It is recommended that each laboratory should establish its own reference interval.

**Quality control**

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

**Performance characteristics****1. Measuring range**

0.05g/L-40g/L

Samples above this concentration should be diluted 1+1 with 0.9% NaCl solution and the result multiplied by 2.

**2. Sensitivity**

When a control serum (14 g/L) is used as a sample, the range of absorbance change should not be smaller than 0.2500.

**Interference**

When Hemoglobin  $\leq 5$ g/L, ascorbic acid (vitamin C)  $\leq 1000$ mg/L, triglyceride  $\leq 10$  mmol/L, heparin lithium  $\leq 200$ U/mL, heparin sodium  $\leq 200$ U/mL, rheumatoid factor  $\leq 800$ IU/mL, no effects on the experimental results are to be found. In order not to affect the test results, please pay special attention to the operation.

**Automation**

Special adaptations for automatic analyzers can be made on request.

**Precautions and warnings****General precautions**

For in vitro diagnostic use only.

Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines.

**Precautions for measurement**

Specimens should be treated as potentially infectious (HIV, Hepatitis B virus, Hepatitis C virus, etc.) and handled with appropriate caution. Reagents with different lot numbers should not be interchanged or mixed.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and

other findings.

Material safety data sheet available for professional user on request.

**Dubai - Global Office****Xportin Global FZE**

Business Centre,

Sharjah Publishing City Free Zone,

Sharjah, United Arab Emirates

Formation No. 4307793.01

Mobile: +971 506938802 / +91 98241 88802

Email: [xportinglobalfzelle@gmail.com](mailto:xportinglobalfzelle@gmail.com)

Website: [www.xportin.com](http://www.xportin.com)

**India - Head Office:****Devasvi Enterprise LLP (Manufacturer)**

Navsari, Gujarat – 396445

Email: [devashvienterprisellp@gmail.com](mailto:devashvienterprisellp@gmail.com)